

PERMEABILITY TEST SET

BS 1377, ASTM 46-47

The Combination Permeameter is designed for laboratory determinations of permeability of either fine-grained or coarse-grained soils. Generally, soils containing 10 percent or more particles passing a no. 200 sieve are tested using the falling head assembly. More granular soils, containing 90 percent or more particles retained on the no. 200 sieve, are tested using the constant head assembly.



SPECIFICATIONS:

Product Code : BSO-215/1 Constant Head Permeability Test Set

Including : 1 m stand with 3 tubes and meter scale, water tank, transparent cells (1 piece 75 mm, 1 piece 114 mm)

Product Code : BSO-215/2 Falling Head Permeability Test Set

Including : 1.5 m stand with 3 tubes and meter scale, water tank, steel cell and soaking tank