

BESMAK – BCE-T/1.0290 DIGITAL BLAIN PERMEABILITY DEVICE



General Specifications:

Electrical blain permeability device, Semi-Automatic

With electric pump and time registration,

Complete with: EN 196-6, BS 4550, ASTM C204

Filter papers 12.8 mm dia. (1000 pcs.)

Fill oil 50 ml, light grease, 50 ml, plug,

Thermometer, brush and funnel

Dim.(wxdxh) = 170 x 300 x 410 mm

Weight: 4.6 kg