

CONE PENETROMETER FOR LIQUID LIMIT

BS 1377 : 2

The cone penetrometer test method for Liquid Limit is based on the relationship between the moisture content and the penetration of a cone into a soil sample. This method eliminates test results dependent upon operator skills and provides a visual measurement of penetration.

Comes with needles which are according to specifications.

SPECIFICATIONS:

Product Code	: BSO-202/2
Weight	: 17 kg
Dimensions	: 250 x 380 x 574 mm



PLASTIC LIMIT TEST SET

ASTM D558 - ASTM D559 - ASTM D698 - ASTM D1557 - ASTM D1883

The plastic limit determines the lowest moisture content of a soil, by which a sample can be rolled into threads 3 mm. dia. without breaking the same neither longitudinally or transversely.

SPECIFICATIONS:

Product Code	: BSO-203
--------------	-----------

Accessories

- 1 piece of 300 x 300 mm glass plate
- 1 piece of 3 mm diameter x 20 cm length stainless steel rod
- 1 piece of Porcelain pot
- 1 piece of elastic spatula
- 6 pieces of 75 mm diameter x 30 mm height sample box

