

BULK DENSITY

EN 1097-3, EN 12350-6, ASTM C138 C29, BS 812

Used to determine the bulk density of aggregates.

Coated with fully electrostatic dye. Produced different volumes according to the diameter x height ratio given in standard.

Capacities are 1-3-5-7-10-15-20 and 30 lt.



SPECIFICATIONS:

Product Code	:	BGE-622/1 1 lt.	BGE-622/5 10 lt
		BGE-622/2 3 lt	BGE-622/6 15 lt
		BGE-622/3 5 lt	BGE-622/7 20 lt
		BGE-622/4 7 lt	BGE-622/8 30 lt


	Weight	Width (mm)	Length (mm)	Height (mm)
1 lt	: 2 kg	100	150	130
3 lt	: 4 kg	165	250	165
5 lt	: 4,5 kg	165	250	270
7 lt	: 5 kg	221	310	206
10 lt	: 8 kg	221	310	295
15 lt	: 12 kg	265	355	300
20 lt	: 12 kg	289	379	335
30 lt	: 18 kg	370	460	300

PARTICLE DENSITY AND WATER ABSORPTION (ABRAHAM CONE AND TAMPER) TEST SET

EN 1097-6, BS 812, ASTM C198, AASHTO T84

Used to determining the specific gravity and absorption of fine aggregates smaller than 2 mm. Apparatus is manufactured from plated steel for protection against corrosion.

Consist of Abraham Cone and tamper.



SPECIFICATIONS:

Product Code	:	BAG-413
Weight	:	500 gr