

MELTING POT

It is used to melt the capping compound, paraffin and other materials.

SPECIFICATIONS:

Product Code	: BCO-132/1 Besmak	BCO-132/2 Impact	BCO-132/3 Humboldt
Capacity	: 3 L	6 L	3,8 L
Temperature Range	: 0-250 °C		
Power	: 220 V / 50 Hz		
Weight	: 9 kg		



CAPPING OF CYLINDERS

EN 12390-3, ASTM C31, C192, C617, AASHTO T23, T126

It is used for flattening cylindrical sample surfaces.

Comes with standard apparatus for 75, 100, 150 mm diameter.

SPECIFICATIONS:

Product Code	: BCO-118
Weight	: 12 kg
Dimensions	:



ELECTRICAL CORE DRILLING MACHINE

EN 12404-1

This is a method of assessing hardened construction materials such as concrete is by cutting a specimen from the area to be checked.

This core drilling equipment is designed to cut cores up to 150 mm diameter from concrete and similar hard construction materials.

The equipment comprises a base unit, support column, feed unit and electrical motor.

It has variable size and power models and comes with stand, water spray apparatus and necessary tools.

Available in two models such as 360° tilting or fixed.

Core drill bits are optional!

SPECIFICATIONS:

Product Code	: BCO-117/1	Electrical Model 360° tilting ,Milwaukee motor
	BCO-117/101	Stand
	BCO-117/2	Electrical Model fixed, Milwaukee motor
	BCO-117/3	Electrical Model fixed, Bosch motor
	BCO-117/301	Stand
Bits	: BCO-117/4	Drill Bit Ø 75 mm
	BCO-117/5	Drill Bit Ø 100 mm
	BCO-117/6	Drill Bit Ø 150 mm
Core Drill Height:	10-50 mm	
Speed	: Two speed 450 rpm , 900 rpm	
Power Supply	: 220 V AC, 50 Hz	
Weight	: 85 kg	
Dimensions	: 700 x 410 x 1200 mm	

